

MICROGEN BIOPRODUCTS LTD

FrozenFresh Cells™ Cell Lines for Viral Culture

Mixed and Single Cell Lines in Frozen Ampoules

FrozenFreshCell R-Mix ampoule (A549+Mv1Lu)

- To seed 50 shell vials only
- Rapid results for DFA detection of respiratory viruses
- Supports growth of respiratory and some enteric viruses
- Available as FreshCell and ReadyCell shell vials

FrozenFreshCell R-Mix *Too* ampoule (A549+MDCK)-to seed 50 tubes

- Excellent replacement for Primary Rhesus Monkey Kidney Cells
- Excellent CPE and virus production at least equivocal PMK's
- Good results for DFA detection of respiratory viruses (48 hours)
- Supports growth of respiratory viruses and some enteric viruses

FrozenFreshCell *Super* E-Mix ampoule (A549+BGMK+hDAF)- to seed 50 tubes

- Optimal mixed cell combination for enteric viruses detection and isolation
- Can be used in conjunction with R-Mix or R-Mix *Too* to cover respiratory and enteric viruses

FrozenFreshCell MRC-5 ampoule - Human Embryonic Lung Fibroblasts—to seed 50 tubes

- Can be used in conjunction with R-Mix or R-Mix *Too* to detect and isolate respiratory and enteric viruses

FrozenFreshCell Hep-2 ampoule –Human epidermoid carcinoma—to seed 50 tubes

- Indicated for influenza detection and isolation
- Can be used in conjunction with other single cell lines to detect and isolate respiratory viruses detection and isolation

FrozenFreshCells are supplied in cryovials containing sufficient cells for seeding 50 13mm tubes or shell vials (R-Mix only).

The Cryovials can be stored at –70°C for 6 months or in liquid nitrogen indefinitely.

Protecting Food and Health


MICROGEN BIOPRODUCTS LTD

1 Admiralty Way
Camberley
Surrey
GU15 3DT
United Kingdom


Phone: +44 (0)1276 600081
Fax: +44 (0)1276 600151
E-mail: sales@microgenbioproducts.com
www.microgenbioproducts.com



Shell Vials—R-Mix



Tubes R-Mix Too, Super E-Mix, Single Cells



Part Codes:

96-00050 R-Mix Cryovial – to seed 50 shell vials with coverslips

97-00050 R-Mix *Too* cryovial – to seed 50 13mm tubes

92-00050 *SuperE*-Mix cryovial – to seed 50 13mm tubes

51-00050 MRC-5 cryovial – to seed 50 13mm tubes

83-00050 MDCK cryovial – to seed 50 13mm tubes

57-00050 Hep-2 cryovial – to seed 50 13 mm tubes

84-00100 Vero Cells

10-200100 Planting medium

10-330100 R-Mix Refeed Medium

10-380100 *SuperE*-Mix Refeed Medium

10-320100 Single Cell Refeed Medium