

MICROGEN BIOPRODUCTS LTD



UNIVERSAL TRANSPORT MEDIUM (UTM)

SEE OTHER SIDE FOR FLOCKED SWABS

Convenient Transport Solution

Room temperature storage for ease of use

Maintain Multiple organisms in one media

UTM has been shown to transport viruses, Chlamydia, Mycoplasma and Ureaplasma exceedingly well.

Flexible packaging, multiple fill volumes

Select from multiple formats of patented Flocked Swabs, packaged with 1mL or 3mL volumes of UTM


Protecting Food and Health

MICROGEN BIOPRODUCTS LTD

1 Admiralty Way
Camberley
Surrey
GU15 3DT
United Kingdom

Phone: +44 (0)1276 600081
Fax: +44 (0)1276 600151
E-mail: sales@microgenbioproducts.com
www.microgenbioproducts.com





FLOCKED SWABS

IMPROVE ASSAY SENSITIVITY

FOR ALL SPECIMEN COLLECTION NEEDS

Increased Sensitivity

Patented flocking technique increases surface area of the swab and puts the collection material where it needs to be.

Decrease rejected samples

Use of flocked swabs increases the quality of samples taken, whilst reducing the number of unacceptable samples.

Flexible packaging, multiple fill volumes

Available in multiple formats with industry leading Universal Transport Medium.

Available in three formats;
Urethral, Regular and Nasopharyngeal.

| Product Code | Description |
|--------------|---|
| 502CS01 | Regular size swab, plastic applicator, nylon fibre tip, sterile in peel pouch |
| 501CS01 | Urethral minitip, plastic applicator, mylon fibre tip, sterile in peel pouch |
| 503CS01 | Flexible Nasopharyngeal minitip, plastic applicator, nylon fibre tip, sterile in peel pouch. |
| 330CHL | 3 mL UTM-RT medium in 10mL size screw cap with internal shaped conical bottom. Each tube contains three-3 mm glass beads. |
| 350CHL | 1 mL UTM-RT medium in 10mL size screw cap with internal shaped conical bottom. Each tube contains three-3 mm glass beads. |
| 401C (404C) | 3 mL (1 mL) UTM with Urethral minitip plastic applicator nylon fibre tip, sterile in peel pouch |
| 402C (405C) | 3 mL (1 mL)UTM with Regular size swab, plastic applicator nylon fibre tip, sterile in peel pouch |
| 403C (406C) | 3 mL (1 mL) UTM with flexible Nasopharyngeal minitip plastic applicator, nylon fibre tip, sterile in peel pouch. |